



		エッペンドルフ												ギルソン																	
		リファレンス4910						リサーチV						リサーチM			ピペットマンP					マルチ									
容量	品番	0.1~2.5μL	0.5~10μL	2~20μL<グレー・ボタン>	2~20μL<イエロー・ボタン>	10~100μL	50~200μL	50~250μL	100~1000μL	0.1~2.5μL	0.5~10μL	2~20μL	10~100μL	20~200μL	100~1000μL	500~5000μL	1000~10000μL	0.5~10μL	10~100μL	30~300μL	P-2/0.1~2.0μL	P-10/0.5~10μL	P-20/2~20μL	P-100/20~100μL	P-200/50~200μL	P-1000/200~1000μL	P-5000/1000~5000μL	P-10ml/1000~10000μL	8×20/12×20	8×300/12×300	
0.1~10μL	204CS	●								●								●			●	●									
	207CS	●	●							●								●				●	●								
	10S	●								●												●	●								
1~20μL	20S			●							●												●								
1~100μL	100S				●							●							●					●							
1~200μL	703CS					●							●												●						
	200S					●							●												●						
	801CS					●							●												●						
1~1000μL	1000S							●						●												●					
	1200S							●						●												●					
1000~5000μL	5000S																										●				
	405CS														●																
	407CS																										●				
1000~10000μL	409CS																											●			

