

ハイパーフィルターチップ適合表

●=適合

◎=容量に注意

※この適合表をご参照のうえサンプルにてご確認ください。

		ワトソン																																		
		NEXTYシングル										NEXTYマルチ	NEXTY容量固定						NEXTYオートクレーバブル						NEXTY-S											
容量	品番	NT-2/0.2~2μL	NT-10/1~10μL	NT-20/2~20μL	NT-100/10~100μL	NT-200/20~200μL	NT-1000/100~1000μL	NT-5000/500~5000μL	NTLG-200/20~200μL	NTIG-1000/100~1000μL	NT-810/1~10μL	NT-8200/20~200μL	NT-F10/10μL	NT-F20/20μL	NT-F100/100μL	NT-F200/200μL	NT-F1000/1000μL	NT-F5000/5000μL	NT-A2/0.2~2μL	NT-A10/1~10μL	NT-A20/2~20μL	NT-A100/10~100μL	NT-A200/20~200μL	NT-A1000/100~1000μL	NT-A5000/500~5000μL	NT-S2/0.2~2μL	NT-S10/1~10μL	NT-S20/2~20μL	NT-S100/10~100μL	NT-S200/20~200μL	NT-S1000/100~1000μL	NT-S5000/500~5000μL				
0.1~10μL	125-10S	●	●							●		●							●	●							●	●								
	125-10XS	●	●							●		●								●	●							●	●							
1~20μL	125-20S			●										●								●								●						
1~100μL	125-100S				●										●								●							●						
1~200μL	125-200S					●			●		●					●								●							●					
	125-200XS					●			●		●					●								●							●					
1~1000μL	125-1000S						●		●								●							●							●				●	
	125-1200S						●		●								●							●							●				●	

		エッペンドルフ														ギルソン															
		リファレンス4910							リサーチV							リサーチM			ピペットマンP					マルチ							
容量	品番	0.1~2.5μL	0.5~10μL	2~20μL<グレー・ボタン>	2~20μL<イエロー・ボタン>	10~100μL	50~200μL	50~250μL	100~1000μL	0.1~2.5μL	0.5~10μL	2~20μL	10~100μL	20~200μL	100~1000μL	500~5000μL	1000~10000μL	0.5~10μL	10~100μL	30~300μL	P-2/0.1~2.0μL	P-10/0.5~10μL	P-20/2~20μL	P-100/20~100μL	P-200/50~200μL	P-1000/200~1000μL	P-5000/1000~5000μL	P-10ml/1000~10000μL	8×20/12×20	8×300/12×300	
0.1~10μL	125-10S		●								●							●			●	●									
	125-10XS		●	●						●								●			●	●									
1~20μL	125-20S				●						●												●								
1~100μL	125-100S					●						●						●						●							
1~200μL	125-200S						●						●												●						
	125-200XS						●						●												●						
1~1000μL	125-1000S							●						●												●					
	125-1200S							●						●												●					

